

Air Mattress Cell

Model: AM40





- Cell mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Comfortable Nylon mattress
- High quality pump
- Weight indicator pressure adjustment setting
- 12 minutes cycling time
- Near silent low noise
- For treatment of patients at medium to high risk of pressure ulcers



Adjustable hangers to suit any kind of hospital bed frame


2 to 1 alternating therapy




Low noise design



Comfortable Nylon mattress



2 to 1 alternating therapy



Easy to use adjustment



Adjustable hangers to suit any kind of hospital bed frame

Accessories



High Quality Rubber Tube



fuse



Spare Cell

Model	Qty per carton	Carton volume
AM 40	2 pcs	0.0511cbm/ctn

Storage and Transportation Condition
-18°C to 43°C (0°F to 110°F) 10~90% RH

Distributed by:

CONSTAR X PTY LTD

Unit 4, 68 Parramatta Road, Underwood, QLD 4119, Australia. Ph: 1300 778 924

www.constar.com.au